


bkFusion TR Steps Quick Guide

Keyboard functions with bkFusion for bk3000/bk5000


Patient, Transducer, Preset Modifications

The top controls for altering exam information is inactive when using bkFusion. You will not be able to modify the following:

- ❌ Enter Patient ID
- ❌ Change Transducer
- ❌ Change Preset



If you have any issues or problems using bkFusion, please contact your BK Medical representative.



Biopsy Steps for TR bkFusion on bk3000/bk5000





1 Click on the **Fusion** workflow tab
(bottom right corner of the monitor at top of the workspace)

2 Click **Start Procedure**


3 Click **TRBx**, and then enter the first two letters of the patient's name. A drop-down list under Patient Name will allow you to select the specific patient, and the remaining patient details will auto-populate. Click **Start Procedure** to begin the fusion exam.

4 The most recent **BXplan** will load on top of the ultrasound image in the upper screen (mid-gland contour). If contours are not displayed, click the **Load Contours** button

5 Alignment

- Make sure signal strength is between 0 and 1 for optimal alignment. Move the EM transmitter or readjust patient position to acquire the best strength.
- Switch to the transverse plane to start alignment (press the bottom button on the transducer or  on the keyboard. Move the transducer along its axis until the MRI contour matches the outer boundary of the prostate, Use the   buttons on the right side of the screen and hold the **Select** key to realign the contours. When the ultrasound image and the prostate contour are correctly aligned, click the **Lock Fusion** button on the left side of the screen.
- Switch to the sagittal plane (use the top button on the transducer or  on the keyboard) and repeat the steps above.

6 Click **Next** to continue to the biopsy screen

7 You take the biopsies in the sagittal plane. Press  on the keyboard to enable the biopsy line. You are now ready to biopsy

8 After each biopsy, use the trackball and the **Select** key to place a biopsy marker on core locations

NOTE: To remove a biopsy market, move the trackball cursor precisely on top of the core, press  and click **Remove Core**, or click the **Remove Last Core** button.

9 Once all biopsy specimens are acquired, click the Save button in the Biopsy Review box to save all data. Select the chosen destination and enter a description for the series

10 Click **Display Report** in the Biopsy Review box on the left-hand side of the screen. Select your Structured Reporting Template from the Notifications pop-up window

NOTE: You can save the report to USB: Click **Create PDF**, select destination (USB), enter name for the report file, and click OK.

11 Click **Sessions** and **End All Sessions** to end exam

Configuring bkFusion Hardware Quick Guide

bkFusion is compatible with the USB foot switch, which can be set as 'mouse left-click'. See the bk3000 & bk5000 Advanced User Guide for foot switch configuration information.

Connect the EM Control Unit to the bk3000/bk5000 system using the USB cable (USB ports are located in the rear of the scanning engine). Connect the power cable, then switch the unit ON



Connect the EM Transmitter to the EM Control Unit: Insert the Transmitter cable into the Transmitter socket until an audible 'click' is heard



Connect the Sensor to the EM Control Unit: Insert the Sensor cable into socket 1 until an audible 'click' is heard



E14C4t (9018)


Push the Sensor into the Sensor Clamp channel. The small indentation on the tip of the Sensor must be facing the groove in the Sensor Clamp channel, as indicated



Slide the other Sensor Clamp channel over the steel stud on the side of the transducer until a 'click' is heard




The numbered and colored scale under the ultrasound image on the system monitor represents the Sensor's proximity to the EM transmitter. Green represents the optimal proximity. To mitigate EM interference, the signal should be either 0 or 1




E13C2 (9029)


Open the Sensor Clamp and push the Sensor into the channel, as indicated. The small indentation on the tip of the Sensor must be facing the groove in the Sensor Clamp channel




Align the hole in the Sensor Clamp with the transducer button



Press the two ends of the Sensor Clamp together until a 'click' is heard




Setting up the EM Transmitter Quick Guide




The EM Transmitter box must be positioned at least 6 inches / 15 centimeters from the **sensor on the transducer**, and no further away than 12 inches / 30 centimeters. The "**Eiffel Tower**" marking on the EM Transmitter must point directly at the **sensor on the transducer**. There must be no metal within the 18-inch / 45-centimeter range.

MRI data import from USB/CD Quick Guide



MRI data import Quick Guide



Connect transducer. Select **Prostate** Exam Type

Click the **Fusion** tab


Click the **Advanced View** tab

Requirements for MIMcloud


- A MIMcloud group has been configured for your site by MIM Software
- Users have MIMcloud accounts
- The BK Ultrasound system has Internet access

...from Remote Patient List

Click on the appropriate patient from the Remote Patient List.




Click the **Transfer Patient Studies** tab, then double-click **Primary Patient List**. Click Yes when prompted by the pop-up window



...from MIM-Cloud (US Only)

Double-click the **MIM Cloud** on the left-hand side of the screen



Enter your login details in the pop-up window

Click on the appropriate patient in the MIM-cloud database

Click the **Transfer Patient Studies** tab, then double-click **Primary Patient List**. Click Yes when prompted by the pop-up window